

LongAxis Manual

By Keith R. Carney, Kwan Lab, University of Utah

August 2019




Topics index	Page
A. Starting LongAxis.....	2
B. Loading and processing confocal z-stack for the first time.....	2
C. Select 3D rendering volume.....	2
D. Cell selection.....	3
E. Adding ellipsoid fit and displaying LongAxis vectors.....	3
F. Exporting data and cell index import/export.....	4
G. Specifying reference point and overlay data slices.....	4
H. Creating urchin plots.....	5
I. Filtering.....	5

A. Starting LongAxis

- Open MATLAB.
- Type "longaxis" at the prompt and LongAxis should open.


B. Loading and Processing Confocal Z-stack for the First Time

1. Click on **Navigation** tab.
2. Click on **Process Image Data** button, and **Select Image Stack** window will open.
3. Select file location of 8bit grayscale z-stack.
4. Input pixel size and z-step size in microns and hit <enter>.
5. Click **Process Image Stack** button after it becomes enabled.
6. LongAxis will now start the 2D image processing. *This could take 5-30 minutes depending on data file size and your computer's processing speed.*
7. When 2D processing is done the **Save Data** button will be enabled. Click **Save Data**, and save processed image data. For subsequent analysis of the same data set, you can skip the 2D processing and quickly load the saved processed data set into LongAxis using **Load Data**.





C. Select 3D Rendering Volume

9. Move **Z-slider** until a representative sample of the image is displayed.
10. Click on the **Draw Select** button.
11. Move cursor to data-viewing window and define the XY region of interest by drawing a rectangle (left-mouse click when done).
12. Move Z-slider up and down and set upper and lower z boundaries by pressing **z1** and **z2** buttons (press again to reset).
13. After XY and Z boundaries are set, press the **Render 3D Cells** button.
14. After rendering is completed, click on **Short Cuts** button for a list of key-presses that will make interacting with the rendering faster and easier.





D. Cell Selection

- Clicking on a cell in the rendered object list (1a), or a cell in the rendering figure (1b), will highlight the cell in both locations as well as display it in the cell-display-figure (1d). The blue tinted sphere around the rendered cell represents the Minimum Bounding Sphere Diameter (MBSD) (1c). Pressing <V> on the keyboard will control the visibility of the cell, and <H> controls highlighting. **Visibility defines which cells are included in data exports, histograms, and convergence point calculations.** Currently selected cell states are also displayed (1e).




E. Adding Ellipsoid Fit and Displaying Longaxis Vectors

- Turn on **ellipsoid fit** by selecting "Ellipsoid Fit" under fit selection menu.
- Display ellipsoid **long-axis vectors** by selecting "cylinder".
- When "Ellipsoid Fit" is turned on the cell-display-figure shows: Ellipsoid overlay (**red**), ellipsoid longaxis (**blue line**), and Minimum Bounding Circle Diameter (MBCD) (**black**). The MBCD is calculated in the plane (**pink octagon**) that runs through the cell's center of mass and is perpendicular to the ellipsoid-fit's long axis.
- Ellipsoid-fit stats are shown in blue.
- The zoom control for the cell-display-figure is expressed in terms of the view angle (degrees). *Small angles are far, and large angles are close.*




F. Exporting Data and Cell Index Import/Export

- Pressing **Export Data** will create a text file (easily imported into MS Excel) that contains the rendered cell measurements (**D1e**). If **Ellipsoid Fit (E1)** and **Cylinder (E2)** are turned on, the ellipsoid measurements are included in the exported data file (**E4**). (*note: orientation angles will only be calculated if there is a reference point specified [see G3]*)
- The **Cell Index Export** button will provide a list of all the current visible/selected cells. *It is recommended copying these index numbers into a MS Word document and labeling them according to their region (example: lens, retina,...).*
- The **Cell Index Import** button allows you to input cell index numbers and control visibility options for them.





G. Specifying Reference Point and Overlay Data Slices

- Select **IMAGE** panel.
- Turning on X,Y, and Z slices (**2a**) will create image-data cross-sections overlaid with the rendered cells. Cross-sections also have transparency (**2b**), smoothing (**2c**), saturation (**2d**), and colormap (**2e**) controls.
- A reference point for orientation angle calculation can be calculated in three different ways:
 - If **Ellipsoid Fit (E1)** and **Cylinder (E2)** are turned on, the **Calculate Convergence Point (3a)** button is enabled. The calculation uses all of the current visible cells.
 - Pressing **Record Current Point (3b)** will record the point of intersection of all three image-data slices.
 - Lastly, the reference point can also be input manually (**3c**).



H. Creating Urchin Plots

1. Select **IMAGE** panel.
2. Turn on Ellipsoid Fit (**E1**) and Vectors (**E2**) and specify a reference point (**3a, 3b, or 3c**).
3. Press **Build Urchin Plot** button.



I. Filtering

1. The visibility of cells can be filtered based on several parameters. *note: all filter functions (except MBSD) require the Ellipsoid Fit (E1) and longaxis vectors (E2) to be turned on.*
 - a. **MBSD**: Filter cells based on their Minimum Bounding Sphere Diameter.
 - b. **r1/r2** ratio: cells can be filtered by the ratio of the fitted ellipsoid's r1 and r2 axis length.
 - c. **r2/r3** ratio: (Same as r1/r2 except the other axes combination)
 - d. **PVD**: Filter cells based on the percent volume difference of the cell and fitted ellipsoid.
 - e. **MBSD/MBSD** ratio: Filter cells by the ratio of their Minimum Bounding Sphere Diameter and Minimum Bounding Circle Diameter.

